



Pressure Care Birthing Mattresses

Made to order

STM Healthcare birthing mattresses are designed to incorporate all of the benefits of pressure care mattresses plus the extra security of welded seams and protected zips to prevent the ingress of fluids. Mattresses can be customised to suit most delivery bed frames.



Birthing		Birthing Set HR Welded STM/30B
Ordering code		
Risk category	See note 1	low/medium
Construction		
Sizes		made to order
Layered foams		
Cover fabric		
HR hydrolysis resistant fabric		
Water-resistant zip		
High frequency welded seams		
Anti-slip base		
Cover colour		plum
Weight guide		
General use, see note 2		233 kg 37 stone
Warranty period, years, see note 3		
Foam		1
Cover		1
Optional Constructions		
Profiling interior hinges		
	standard	
	optional	
HR	Hydrolysis resistant fabric coating	
Note 1	Indication of patient risk category - not a prescriptive criteria.	
Note 2	Weight limit can vary depending on type and length of use - these weight limits are to be used as a guide only and in conjunction with clinical judgements. For constant long term use please contact STM Healthcare for information.	
Note 3	Refer to written warranty for full details	
Note 4	Fire retardancy BS7177:2008 medium hazard: Fillings comply with Furniture & Furnishings (Fire) (Safety) Regulations 1988.	
Note 5	Fabrics are multi-stretch, wipe clean, waterproof, vapour permeable, PU coated knitted nylon. Cover colours may vary.	
Note 6	Produced from superior grade CMHR foam grades.	

STM
HEALTHCARE

A DIVISION OF RECTICEL LIMITED 

STM Healthcare

Azalea Close, Clover Nook Industrial Park
Alfreton, Derbyshire DE55 4QX

Tel: 01773 830426 Fax: 01773 830427

stmhealthcare@recticel.com

www.stm-healthcare.co.uk