

Pressure Care Operating Table Top Sets



Operating table top sets have been suspected of contributing to the onset of pressure sores. Contributing factors such as reduced patient temperature leading to reduced circulation, immobility and long periods of time spent on firm surfaces can affect tissue viability. STM Healthcare can provide bespoke table top sets with pressure reducing surfaces and welded seams to help prevent tissue damage.



Operating			
		Std Welded	HR Welded
Ordering code		STM30	STM30E
Risk category	See note 1	low/ medium	low/ medium
Construction			
Size		made to order	
Layered foams		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Non-turning		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Cover fabric			
HR hydrolysis resistant fabric		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
High frequency welded seams		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Cover colour		blue/black	spearmint
Weight guide			
General use, see note 2		185 kg 29 stone	185 kg 29 stone
Warranty period, years, see note 3			
Foam		1	1
Cover		1	1
Optional Constructions			
Profiling interior hinges		<input type="checkbox"/>	<input type="checkbox"/>

standard
 optional

- HR Hydrolysis resistant fabric coating
- Note 1 Indication of patient risk category - not a prescriptive criteria.
- Note 2 Weight limit can vary depending on type and length of use - these weight limits are to be used as a guide only and in conjunction with clinical judgements. For constant long term use please contact STM Healthcare for information.
- Note 3 Refer to written warranty for full details
- Note 4 Fire retardancy BS7177:2008 medium hazard: Fillings comply with Furniture & Furnishings (Fire) (Safety) Regulations 1988.
- Note 5 Fabrics are multi-stretch, wipe clean, waterproof, vapour permeable, PU coated knitted nylon. Cover colours may vary.
- Note 6 Produced from superior grade CMHR foam grades.



A DIVISION OF RECTICEL LIMITED

STM Healthcare

Azalea Close, Clover Nook Industrial Park
Alfreton, Derbyshire DE55 4QX

Tel: 01773 830426 Fax: 01773 830427

stmhealthcare@recticel.com

www.stm-healthcare.co.uk