

Pressure Care Cot & Paediatric Mattresses

Made to order

In certain circumstances children and infants can be at risk of tissue damage. STM Healthcare cot and paediatric mattresses are designed and manufactured using the same techniques and expertise incorporated into pressure care mattresses for adults but using grades of foam suitable for lower body weights. STM Healthcare cot and paediatric mattresses are available in all sizes.



Paediatric				
		Cot Sewn	Crib Sewn	Junior Sewn
Ordering code		STM/LNCot	STM/LNCrib	STM/LNJR
Risk category	See note 1	low/medium	low/medium	low/medium
Construction				
Sizes		made to order		
Layered foams				■
4-way turning		■	■	■
Cover fabric				
Water-resistant zip		■	■	■
Sewn cover seams		■	■	■
Cover colour		spearmint	spearmint	spearmint
Weight guide				
General use, see note 2		38 kg 6 stone	19 kg 3 stone	185 kg 29 stone
Warranty period, years, see note 3				
Foam		8	8	8
Cover		4	4	4

■ standard
□ optional

- HR Hydrolysis resistant fabric coating
- Note 1 Indication of patient risk category - not a prescriptive criteria.
- Note 2 Weight limit can vary depending on type and length of use - these weight limits are to be used as a guide only and in conjunction with clinical judgements. For constant long term use please contact STM Healthcare for information.
- Note 3 Refer to written warranty for full details
- Note 4 Fire retardancy BS7177:2008 medium hazard: Fillings comply with Furniture & Furnishings (Fire) (Safety) Regulations 1988.
- Note 5 Fabrics are multi-stretch, wipe clean, waterproof, vapour permeable, PU coated knitted nylon. Cover colours may vary.
- Note 6 Produced from superior grade CMHR foam grades.

STM
HEALTHCARE

A DIVISION OF RECTICEL LIMITED 

STM Healthcare
Azalea Close, Clover Nook Industrial Park
Alfreton, Derbyshire DE55 4QX
Tel: 01773 830426 Fax: 01773 830427
stmhealthcare@recticel.com
www.stm-healthcare.co.uk