



'Very high dependency' mattresses

5 years continuous clinical trials

The Prima mattress is manufactured using soft, high density foam layers for outstanding comfort, support and durability. It provides exceptional pressure area care and is available for all types and sizes of bed base.

In response to the need for reduced levels of lifting and handling, non-turning versions of the Prima mattress are available, one of which includes a superior high density visco elastic foam insert.

- Suitable for profiling and non-profiling bed frames
- Unique machine-hinged profile design (optional)
- Does not impede CPR
- 8 year warranty for the foam*
- 4 year warranty for the cover*
- Standard size 193 x 91 x 15 cm - alternative sizes available on request (welded and sewn seams)
- Fully covered zip fastening on non-turn models
- Vapour permeable, multi-stretch fabric cover
- Anti-slip base fabric on non-turn models
- White substrate fabric to help identify interior staining
- Complies with BS 7177:2008 medium hazard
- Audit box, rotation cycle and care and cleaning instructions printed on to the cover
- Weight guide - up to 37 stone/233 kg depending on model



* refer to written warranty details

Prima Mattresses

	Std Sewn	HR Sewn	HR Welded	HR Welded Non-turn	HR Welded Non-turn Visco STM/6V
Ordering code	STM/6	STM/6E	STM/6W	STM/6N	STM/6V
Risk category See note 1	very high	very high	very high	very high	very high
Construction					
193 x 91 x 15 cm	■	■	■	■	■
Layered foams	■	■	■	■	■
4-way turning	■	■	■		
Non-turning				■	■
Cover fabric					
HR hydrolysis resistant fabric		■	■	■	■
Water-resistant zip	■	■	■	■	■
Sewn cover seams	■	■			
High frequency welded seams			■	■	■
Anti-slip base				■	■
Cover colour	plum	plum	plum	plum	plum
Weight guide					
General use, see note 2	213 kg 34 stone	213 kg 34 stone	213 kg 34 stone	233 kg 37 stone	171 kg 27 stone
Warranty period, years, see note 3					
Foam	8	8	8	8	4
Cover	3	4	4	4	4

Optional Constructions

Profiling interior hinges

■ standard
□ optional

HR Hydrolysis resistant fabric coating

Note 1 Indication of patient risk category - not a prescriptive criteria.

Note 2 Weight limit can vary depending on type and length of use - these weight limits are to be used as a guide only and in conjunction with clinical judgements. For constant long term use please contact STM Healthcare for information.

Note 3 Refer to written warranty for full details

Note 4 Fire retardancy BS7177:2008 medium hazard: Fillings comply with Furniture & Furnishings (Fire) (Safety) Regulations 1988.

Note 5 Fabrics are multi-stretch, wipe clean, waterproof, vapour permeable, PU coated knitted nylon. Cover colours may vary.

Note 6 Produced from superior grade CMHR foam grades.



A DIVISION OF RECTICEL LIMITED

STM Healthcare

Azalea Close, Clover Nook Industrial Park
Alfreton, Derbyshire DE55 4QX

Tel: 01773 830426 Fax: 01773 830427

stmhealthcare@recticel.com

www.stm-healthcare.co.uk