

Pressure Care Trolley Mattresses



Tissue damage leading to pressure sores is not solely confined to patients in bed or sitting. Many patients spend a great deal of time on trolleys either on entry to the hospital or because they are undergoing lengthy procedures. Not all hospital trolley mattresses offer comfort or pressure care. It makes sense therefore to use a pressure relieving trolley mattress to close this loophole in pressure area care. STM Healthcare trolley mattresses are fully clinically evaluated* and available in all sizes.

* "Comparison of pressure distribution in hospital trolley mattresses", J. Tissue Viability, 1993 Vol 3 No 1.



Trolley

	Std	HR	HR
	Sewn	Sewn	Welded
Ordering code	SMT15	STM15E	STM15EW

Risk category	See note 1	low/medium	low/medium	low/medium
---------------	------------	------------	------------	------------

Construction				
Sizes	made to order			
Profiled surface	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	

Cover fabric				
HR hydrolysis resistant fabric		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Water-resistant zip	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	
Sewn cover seams	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>		
High frequency welded seams				<input checked="" type="checkbox"/>
Cover colour	blue/black	spearmint	spearmint	

Weight guide				
General use, see note 2	140 kg 22 stone	140 kg 22 stone	140 kg 22 stone	

Warranty period, years, see note 3				
Foam	1	1	1	
Cover	1	1	1	

Optional Constructions				
Profiling interior hinges	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

standard
 optional

- HR Hydrolysis resistant fabric coating
- Note 1 Indication of patient risk category - not a prescriptive criteria.
- Note 2 Weight limit can vary depending on type and length of use - these weight limits are to be used as a guide only and in conjunction with clinical judgements. For constant long term use please contact STM Healthcare for information.
- Note 3 Refer to written warranty for full details
- Note 4 Fire retardancy BS7177:2008 medium hazard: Fillings comply with Furniture & Furnishings (Fire) (Safety) Regulations 1988.
- Note 5 Fabrics are multi-stretch, wipe clean, waterproof, vapour permeable, PU coated knitted nylon. Cover colours may vary.
- Note 6 Produced from superior grade CMHR foam grades.



A DIVISION OF RECTICEL LIMITED

STM Healthcare
 Azalea Close, Clover Nook Industrial Park
 Alfreton, Derbyshire DE55 4QX
 Tel: 01773 830426 Fax: 01773 830427
 stmhealthcare@recticel.com
 www.stm-healthcare.co.uk