

Pressure Care Cushions


Neo Range

Neco



An affordable modular foam cushion with profiled upper surface and supportive base.

Neogel



The Neogel cushion uses high resilient foam with a profiled upper surface to increase pressure distribution, encasing a silicon gel sac, which dissipates pressure in areas of high loading.

Neoform




High density foams for greater levels of support means that the Neoform is suitable for heavier patients. Its layered construction provides optimum comfort.

Neotherm



Suitable for heavier patients in either wheelchairs or day care chairs, the combination of high density high resilient foam base and visco elastic heat sensitive foam surface provides support and pressure distribution.

Neotek



A superior performance cushion, the unique 'foam spring' design aids ventilation, durability and mobility. The premium foam grades are ideal for heavier patients.

Neo Cushions

		Neco	Neoform	Neotek	Neogel	Neotherm
Ordering code		STM/N	STM/NF	STM/NK	STM/NG	STM/NT
Risk category	See note 1	low/medium	medium/high	high	high	high
Construction	Safeguard® CMHR foam	●	●	●	●	●
	Fire retardancy, see note 2	●	●	●	●	●
	Incorporates heat sensitive visco elastic foam					●
Cover fabric	Multi stretch, wipe clean, waterproof, vapour permeable, PU coated knitted nylon	●	●	●	●	●
	Welded cover seams	●	●	●	●	●
	Anti slip base	●	●	●	●	●
	Zip protected by flap	●	●	●	●	●
Weight guide	Maximum	95kg 15 stone	127 kg 20 stone	127 kg 20 stone	121 kg 19 stone	127 kg 20 stone
	Colour	blue	blue	blue	blue	blue
Standard sizes	41 x 41 x 8 cm	●	●	●	●	●
	43 x 43 x 8 cm	●	●	●	●	●
	46 x 41 x 8 cm	●	●	●	●	●
	46 x 46 x 8 cm	●	●	●	●	●
	48 x 48 x 8 cm	●	●	●	●	●
	51 x 41 x 8 cm	●	●	●	●	●
Other sizes to special order	51 x 46 x 8 cm	●	●	●	●	●
	FR terry towelling cover	○	○	○	○	○
	Waterproof interliner	○	○	○	○	○
Warranty period	Months	12	12	12	12	12

● standard

Note 1

Note 2

○ optional - further details available on request

Indication of risk - not a prescriptive criteria

BS 6807 : 2006; Fillings comply with: Furniture and Furnishings (Fire) (Safety) Regulations 1988

Neo cushions have been tested at the Department of Orthopaedic Mechanics, University of Salford, Lancs.

STM Healthcare

Azalea Close, Clover Nook Industrial Park
Alfreton, Derbyshire DE55 4QX

Tel: 01773 830426 Fax: 01773 830427

stmhealthcare@recticel.com

www.stm-healthcare.co.uk

STM
HEALTHCARE

A DIVISION OF RECTICEL LIMITED 